# BIOSTIM SPT®

800 990053  
+ 39 059 699600

+ 49 89 59081236

+ 44 1937 547065

+ 39 059 699600

**IGEA**<sup>®</sup>  
CLINICAL BIOPHYSICS

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**Begin using the device only after having carefully read  
this user manual**

## INDICE

1. INTRODUCTION .....	15
1.1 What is BIOSTIM SPT and how it works .....	15
1.2 Therapeutic treatments provided by BIOSTIM SPT .....	15
2. BIOSTIM SPT COMPONENTS .....	15
2.1 Generator description .....	16
2.2 Coil connection to the generator .....	16
3. USE OF BIOSTIM SPT .....	16
3.1 Therapy execution .....	16
3.2 Coil Positioning .....	17
3.3 Recommendations for the correct use of the coil .....	17
3.4 Treatment times .....	18
4. TROUBLE SHOOTING .....	18
5. SAFETY INSTRUCTIONS .....	19
5.1 Recommendations .....	19
5.2 Contraindications and adverse effects .....	20
5.3 Electromagnetic compatibility .....	20
5.4 Biological safety .....	21
6. MANUFACTURER'S RESPONSIBILITY .....	21
7. TECHNICAL DATA .....	22
7.1 Table of Electromagnetic Emission and Immunity Levels .....	23
7.2 Immunity to proximity fields from RF wireless communication equipment .....	24
8. SYMBOLS .....	24

## INTRODUCTION

### 1.1 What is BIOSTIM SPT and how it works

BIOSTIM SPT is an electrical osteogenetic stimulator, which uses electromagnetic fields.

BIOSTIM SPT is a therapeutic appliance and is to be used only on medical prescription.

The system basically consists of an electromagnetic field generator pulsated at low frequency and consists of an impulse signal at a frequency of 75 Hz, with a 10% duty cycle, that is, an activation time of approx. 1.33 milliseconds.

This electromagnetic field has the capacity of inducing an average electrical field of 0.04mV/cm into the bone tissue, which consists in the active component of the signal, able to boost osteoblast activity.

The electrical field is concentrated at the point to be treated by means of coils of suitable shape.

The coil/coils to be applied to the part requiring treatment, provide an extremely homogeneous field and therefore do not require application to be perfectly centred on the part to be treated; this means that the patient can undertake treatment himself, without requiring the supervision of either doctors or nurses.

The generator is controlled by a processor which constantly monitors the efficient function of the appliance, providing rapid warning of any eventual anomalies or malfunctioning during the treatment. For this purpose the appliance is equipped with simple and effective visual and acoustic indicators.

### 1.2 Therapeutic treatments provided by BIOSTIM SPT

BIOSTIM SPT is used principally in orthopaedics and traumatology fields for the treatment of the followings:

- Non-union
- Recent fracture
- Osteotomy
- Necrosis
- Painful Prosthesis
- Risk Fracture
- Osteoporosis Fracture
- Revision Prosthesis
- Bone grafts
- Skin wound healing

## 2. BIOSTIM SPT COMPONENTS

BIOSTIM SPT system consists of: generator (①), coil (②) applied part of the device, power supply (③).

The system is usually supplied with a velcro strip, which can be used to maintain the coil in correct position during treatment. Nevertheless the patient may replace it with other similar types such as elastic straps or bands as required by the user.

## 2.1 Generator description


The generator is complete with:

- A **green** led “” indicating that the appliance is ON and correctly connected to the power supply.
- A **green** led “” indicating magnetic field presence. Under correct conditions of use, the led is on, indicating the generator is providing the right driving signal to the coil.
- A **red** led “” indicating fault/malfunction of the coil; this indicator is accompanied by an acoustic alarm.
- A **yellow** flashing light “” which lights up in the event of the faulty function of the power supply.
- On/off **switch** positioned to the right of the generator (④).
- A socket positioned on the lower side of the appliance “—●—” (⑤) used for connection to the power supply.

## 2.2 Coil connection to the generator

- **METALLIC CONNECTOR:** Connect the coil to the generator by inserting the male circular connector to the female circular connector, which exits from the generator, then turn the two metallic rings present on both the connectors in an opposite sense until the connection is thoroughly tightened.

- **PLASTIC CONNECTOR:** connect the coil connector to the generator connector, push the connector until at "clik" in this way the connection is thoroughly tightened. To separate the coil from the generator, lift the tab on the connector which exits from the generator and exert a slight pull.

## 3. USE OF BIOSTIM SPT

### 3.1 Therapy execution

- Connect the cable of the power supply to the socket on the lower part of the generator “—●—” then connect the power supply plug to a main socket.

- Position the coil on the part to be treated, securing it using the velcro strap provided, as described in **3.2**.
- Switch on the generator positioning the switch to "I" (ON); the device performs an initial test during which all the lights and the acoustic alarm switch on in sequence. At the end of the test, the device starts normal operation.
- **The green light "💡" will come on indicating that the device is correctly supplied; in the same way the 💡 green lights will come on indicating the coil driving signal presence.**
- Undertake the treatment as described in paragraph **3.4**, or as indicated by the doctor who prescribed the treatment.
- Once the application is complete switch off the generator positioning the switch to "O" (off) and disconnect the power supply from the power socket.

### 3.2 Coil Positioning

Position the coil, ensuring that the injury, or the part of the body to be treated, **is at the centre of the area occupied by the coil**; it is recommended positioning the coil on a garment or on a cloth.

Double coils should be placed in opposed way, laterally if possible, so that the treatment area is between the two coils.

Secure the coil using the velcro strap provided or by means of another suitable means (an elastic band for example); **do not tighten or force the contact between the coil and the treatment area.**

Below are some examples of positioning.

Positioning examples with single coil

Positioning examples with double coil

### 3.3 Recommendations for the correct use of the coil

- The coil does not have to be applied directly on the skin; for hygiene purposes, it is **recommended to always place the coil over light clothes**, which cannot affect treatment in any way. In particular if the coil is to be applied on injured skin, wear a light item of clothing that covers the area; if this is not possible, disinfect the coil before applying it to the injured skin.
- For effective treatment is sufficient to fix the coil on the treatment site using the velcro strap provided; **do not force the contact** with the coil.
- Double coils should be placed in opposed way, laterally if possible, so that the treatment area is between the two coils.

- Do not use the coil if the cable connecting it to the generator shows obvious sign of wear or damage.
- The use of the coil under heavy blankets limits heat dispersion and can lead to overheating of the coil surface; in case of excessive overheating, discover the coil. In any case, the generator and the power supply must not be covered during the operation, to allow its aeration.
- Before proceeding with cleaning the coil ensure that the appliance is switched off and disconnected from the power supply. The coil can be cleaned using neutral detergents, never use aggressive detergents or solvents for cleaning it. Regularly disinfect the coil, particularly in the case of prolonged use or direct contact with the skin.
- The wear of the coil coating due to the use does not affect the effectiveness of the therapy. In case of loss of integrity of the coating, the coil must be replaced.

### 3.4 Treatment times

The user should perform the treatment by the number of days specified by the doctor who prescribed the therapy.

BIOSTIM SPT treatment time is of 6/8 hours a day.

Daily therapy should preferably all be in one session; if you prefer to split daily therapy time into several parts, continuous periods should be of **no less than 2 hours**. As there are no adverse side effects, therapy can also be carried out during sleeping hours.

## 4. TROUBLE SHOOTING

BIOSTIM SPT can detect and signal any conditions preventing correct operation. The table describes the system messages and actions to restore correct operation.

IF	THEN
The green light “  ” on the generator fails to come on.	<p>Check if:</p> <ul style="list-style-type: none"> <li>- the generator is connected to the power supply,</li> <li>- the power supply is connected to the main,</li> <li>- the generator is on.</li> </ul>
The red <b>flashing</b> light “  ” comes on, a <b>acoustic flashing signal</b> come out and the green lights  indicating the magnetic field presence is off.	<ul style="list-style-type: none"> <li>- The coil is incorrectly connected to the generator connector.</li> <li>- If, after switching the device off and on again, the coil is correctly connected to the generator, but the failure continues, the coil is faulty and needs to be replaced.</li> <li>- External interference or short interruptions in the mains power supply could cause the suspension of the treatment, signaled by the device as described here on the side. In this case, switch the device off and on again to restore normal operation.</li> </ul>

The green light "  " is on at the same time as the yellow light "  "	- Power supply generates inadequate power and needs to be replaced. - The coil has short-circuited causing a lowering in the power voltage; replace the coil.
Coil replacement	To replace the coil contact the manufacturer or an authorized assistance centres. For the connection of the new coil refer to the instructions in par <b>Errore. L'origine riferimento non è stata trovata.</b>
The red <b>fixed</b> light "  " comes on, a continuous acoustic signal come out and the green lights  indicating the magnetic field presence is off.	The device must be delivered to IGEA assistance centre for maintenance.
<p> <b>The device technical assistance is exclusive pertinency of Manufacturer IGEA S.p.A.</b> In case of failure, or in any case where repairs is necessary, the user should contact a service centre authorized by IGEA S.p.A.</p> <p><b>IGEA S.p.A. - Via Parmenide 10/A - 41012 Carpi (MO) - ITALY</b>  <b>Tel. +39 059 699 600 Fax. + 39 059 695 778 e-mail: <a href="mailto:info@igeamedical.com">info@igeamedical.com</a></b></p>	

## 5. SAFETY INSTRUCTIONS

### 5.1 Recommendations

Closely follow these instructions to ensure optimal and safe operation of the device:

#### Make sure you read the instruction manual before using the device

- BIOSTIM SPT must be used by people who understand and put into effect independently the instructions provided in this manual: otherwise, and if used on children, BIOSTIM SPT can be used under the supervision of people able to understand and put into effect the instructions provided in this manual.
- Keep the device away from children and pets, if any.
- Attention, connecting cables could cause strangulation hazard if incorrectly used.
- Never use the appliance in presence of inflammable gases
- Ensure that the power voltage corresponds to that indicated on the power supply data plate.
- Do not connect any part of unit to other equipment or devices.
- Do not connect to BIOSTIM SPT any part not intended for use and not supplied by the Manufacturer.
- Use only the power supply provided by the manufacturer; the use of different devices could cause damage to the system or in the worst case to the user, for which the manufacturer declines all responsibility.
- Always check that the mains connection cable, cables and connectors are not visibly damaged: in the event of damage contact IGEA or an authorized assistance centre.
- Never handle any of the components of the BIOSTIM SPT system with wet hands in particular during the connection to the mains.
- Never immerse any of the system components in water or liquids of any type and never pour any liquids on to the same.
- Do not cover the generator or the power supply, for example with blankets, during the delivery of therapy to allow their aeration.
- Using the coil under heavy blankets, you may experience a overheating of the coil surface: in this case, perform the therapy without covering the coil.
- It is advisable not to use the system in environments where the temperature is above 30°C, to avoid overheating of the coil surface; if this is not possible, and the temperature of the coil does

cause discomfort, it is recommended to split the daily time therapy in multiple sessions lasting no less than two hours each.

- The device must be used in compliance with the normal rules of hygiene, and should be cleaned regularly. The presence in the environment of use of hair, hair, dust, as well as exposure to direct sunlight, while not causing malfunction of the device, should be avoided.
- Before proceeding with cleaning the coil or the generator ensure that the appliance is switched off and disconnected from the power supply. All BIOSTIM SPT components can be cleaned using neutral detergents; never use any aggressive solvents or detergents for cleaning.
- The coil is single patient
- The Velcro strap can be washed like any garment.
- Avoid knocking the appliance during transportation and movement.
- The device is equipped with mechanisms of self-control of the correct functioning of any anomaly that may occur is indicated by the system and is described in the instruction manual. Any malfunctions therefore do not affect therapeutic efficacy.
- IGEA recommends to repeat safe test at regular intervals not exceeding 24 months of use, to ensure that safety standards are consistently maintained. IGEA, in agreement with the client, is able to perform the recommended safe test.
- IGEA recommended a routine maintenance procedure for the operating parameters check at regular intervals not exceeding 24 months of use to ensure reliable performance; contact IGEA or an authorized assistance centre.

### **5.2 Contraindications and adverse effects**

In the many international studies carried out, no negative side effects have been detected due to the therapy, nevertheless it is necessary to follow the below mentioned precautions:

- In case of pregnancy, certain or supposed, even if no negative effects are described relating to the therapy, as a caution it is suggested to avoid to treat the pelvis bones directly. In any case, always contact the doctor who has prescribed the therapy, he will evaluate case by case the need of prosecuting/interrupting it.
- The pacemakers currently in use today are not affected by electromagnetic fields, however those who are fitted with a pacemaker are advised to verify the features of their own pacemaker with their cardiologist, and to avoid any treatments that involve direct coil application to the chest.
- Fewer than 0,2% of patients complained of burning sensations during treatment. In the event that this occurs it is necessary to reduce the treatment time to 3 hours in the first 7 days of treatment, to be divided into 3 sessions of one hour each; then gradually increase the treatment time by one hour a day until the standard application time is reached. The burning sensation ceases on suspension of treatment.
- The use of the device at the same time the use of drugs not normally presents no contraindication. The doctor who prescribes the therapy should be informed by the patient for any medication.

### **5.3 Electromagnetic compatibility**

BIOSTIM SPT has been tested and certified according to the electromagnetic compatibility standards for medical devices and it is declared suitable for the home environment.

BIOSTIM SPT can be used together with other electrical or electronic devices, providing they also conform to current standards, without causing interference or receiving disturbances. However, the following general requirements must be observed:

- BIOSTIM SPT must not be used adjacent or overlapped with other devices. If adjacent or overlapping use is required, the medical device must be observed to verify normal operation in the configuration in which it is used;
- BIOSTIM SPT requires special precautions regarding electromagnetic compatibility and must be installed and used in compliance with the electromagnetic compatibility information provided below in this manual;
- BIOSTIM SPT must not be used at the same time as other therapies or applications of electro-medical devices that release energy to the patient's body, particularly if they use high-frequency

signals, as these signals may interact with the operation of BIOSTIM SPT and cause unwanted changes in the therapeutic signal;

- The use of accessories and cables other than those specified and provided directly by the BIOSTIM SPT manufacturer, may result in increased emissions or decreased immunity of BIOSTIM SPT and cause improper operation;
- BIOSTIM SPT may be sensitive to electrostatic discharge with a potential > 8kV; the user must take every precaution to avoid discharges to the device that could cause a operational block or a malfunction.
- Electrostatic discharge with a potential of more than 15kV may damage the signal generator until it causes a permanent failure; in this case the generator still maintains the basic safety required and does not constitute any danger for the user or for third parties. **In this case the generator will have to be replaced with an efficient one.**
- Portable and mobile RF communications devices, including peripheral devices such as antenna cables and external antennas, should be kept more than 30 cm away from all BIOSTIM SPT components, including cables. Otherwise the performance of the medical device may deteriorate.

**Device lock:** Electromagnetic interference, such as those due to active cell phones or disturbances to the power system, can interfere with the normal operation of BIOSTIM SPT and cause the device to lock. In the event of a **operational block**, indicated by BIOSTIM SPT as a situation of failure or interruption in the delivery of therapy, the device must be turned off and then on again after a few seconds, using the switch on the side of the device to restore normal operation.

Be sure to remove the source of the disturbance before continuing treatment, if this is not possible, increase the distance to the source of disturbance to at least 1 meter.

**Caution:** Using the device in the immediate vicinity of a shortwave or microwave (in operation) therapy device may cause disturbances in the output signal which, although it does not pose a risk to the patient, could nullify the therapeutic effect; therefore it shall be avoided.

#### 5.4 Biological safety

The safety of the treatment with BIOSTIM SPT has been extensively verified; all tests showed the absence of negative effects of treatment.

### 6. MANUFACTURER'S RESPONSIBILITY

IGEA considers itself responsible for the safety, reliability and performance of the BIOSTIM SPT system **providing that:**

1. Any modifications or repairs are undertaken only by personnel authorized by the same.
2. The appliance is not opened or tampered with in any way by the user or other non-authorized persons.
3. The appliance is used in accordance with the operating instructions of the present manual;
4. The power is supplied using only power supply provided by the manufacturer.
5. The power supply is exclusively used for operating the BIOSTIM SPT system according to the instructions of the present manual.
6. The user regularly perform the safe test and the operating parameters check recommended by IGEA.

IGEA company reserves the right to implement updates that improve the system performances without modifying pictures or instructions of the present manual.

For any further information or revisions that user can contact the local dealer or the manufacturer itself:

**Manufacturer & Maintenance Service Centre:**

**IGEA S.p.A. Via Parmenide 10/A - 41012 Carpi (MO) – ITALY**

Tel. +39 059 699 600 Fax. +39 059 695 778 - [info@igeamedical.com](mailto:info@igeamedical.com)

## 7. TECHNICAL DATA

BIOSTIM SPT complies with MDD 93/42 EEC and subsequent amendments and is  marked under the control of IMQ.

BIOSTIM SPT has an expected lifetime of five years.

### **BIOSTIM SPT Generator**

Model	: BIOSTIM SPT
Power supply	: 12 VDC
Max. absorbed power	: 0,8 A.
Classification according to EN 60601-1	: Class II with Applied Part of Type BF
Classification according to MDD 93/42 EEC	: Ila

**Method of use:** Appliance for continuous use not to be used in the presence of anaesthetic mixtures inflammable on contact with air, oxygen, or nitrous oxide.

Unit with casing having a degree of protection IP10.

### **Characteristics of coil's pilot signal & magnetic field intensity:**

Signal type	: Square wave
Signal frequency	: 75 Hz.
Duty Cycle	: 10% (1,33 msec)
Magnetic field intensity	: 20-30 Gauss (peak value)

### **External power supply**

Model*	MENB1020A1203B01	ME20A1203B01
Brand*	AULT - SL Power	SL Power
Input voltage (wide band)	230VAC (100-240)	230VAC (100-240)
Mains frequency	50-60Hz	50-60Hz
Max. absorbed power	0,100 A	0,200 A
Output voltage	12VDC	12 VDC
Max. output voltage	1,5A	1,5 A
Overload protection	Internal self-resetting protection	Protezione continua interna auto-ripristinante
Insulation category	II	II

\*IGEA reserves the right to provide different models of power supply, tested and approved for the system according to the standards EN60601-1, EN 60601-1-2; use only the power supply provided directly by the manufacturer.

### **Environmental condition of use**

Atmospheric temperature	: 0 - 40°C
Relative humidity	: 30% - 75%
Atmospheric pressure	: 700-1060hPa

### **Environmental transportation and storage conditions**

Atmospheric temperature	: -20 - +70°C
Relative humidity	: 10% - 90%
Atmospheric pressure	: 500 - 1060hPa

**Note:** After removing the device from its packaging it is considered that the conditions of use are also valid for transport and storage.

**Warning: This device is not to be used in areas in which an explosion risk exists.**

**End of life disposal**

BIOSTIM SPT and all its parts are subject to separate collection in the manner established by local authorities.

**7.1 Table of Electromagnetic Emission and Immunity Levels**

BIOSTIM SPT taking into account the risk analysis associated with electromagnetic disturbances, complies with all test levels & complies to EN60601-1-2.

BIOSTIM SPT, can be used in the specified electromagnetic environment. The user must ensure that it is used in an electromagnetic environment with the features described below

Emissions tests	Compliance	Electromagnetic Environment
RF emissions - CISPR 11	Group 1	BIOSTIM SPT, generates RF signals exclusively as a result of the operation of internal electronic circuits. Its RF emissions are very low and hardly cause radio interference in nearby equipment.
RF emissions - CISPR 11 Harmonic Emissions EN 61000-3-2	Class B Class A	BIOSTIM SPT, is suitable for use in all environments, including domestic environments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations /flicker emissions EN 61000-3-3	Compliant	

Immunity Test	EN 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment
Electrostatic discharge (ESD) EN 61000-4-2	± 8 kV contact ± 15 kV air	IEC 60601-1-2 Test level	All environments, including domestic environment <sup>(1)</sup> .
Radiated RF EN 61000-4-3	10 V/m from 80 MHz to 2,7 GHz	IEC 60601-1-2 Test level	Any environment, including domestic environment with portable and mobile RF equipment kept as far as possible from BOSTIM SPT, including connection cables. Minimum distance 30 cm <sup>(2)</sup>
Electrical fast transient/burst EN 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output	IEC 60601-1-2 Test level	All environments, including domestic environment.
Surges EN 61000-4-5	± 1 kV between phases ± 2 kV form phase to ground	IEC 60601-1-2 Test level	All environments, including domestic environment.
Conducted disturbances induced by RF fields EN 61000-4-6	3 Vrms from 150 kHz to 80 MHz 6 V - ISM Frequencies and Radio amateur band	IEC 60601-1-2 Test level	Any environment, including domestic environment with portable and mobile RF equipment kept as far as possible from BOSTIM SPT, including connection cables. Minimum distance 30 cm
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	10 ms – 0% a 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 20 ms – 0% a 0° 500 ms – 70% a 0° 5 s – 0%	IEC 60601-1-2 Test level	All environments, including domestic environment.
Power frequency (50/60 Hz) magnetic fields EN 61000-4-8	3 A/m	IEC 60601-1-2 Test level	All environments, including domestic environment.

<sup>1</sup> Electrostatic discharges of potential > 6kV (in contact) and > 8kV (in air) may cause a blockage of the device or an interruption in the delivery of therapy. Potential 15kV discharges could cause a permanent failure as described in the paragraph

<sup>2</sup> The presence of electromagnetic disturbances at high levels, limited to some intervals of frequency, could cause an interruption in the delivery of the signal, signaled by the device. BOSTIM SPT, considering the risk analysis EMC meets all test levels and complies with the EN 60601-1-2 Standard

## 7.2 Immunity to proximity fields from RF wireless communication equipment

The portable or mobile RF communication devices that may be present in the home, such as wireless telephones, mobile phones, devices for wireless connection to the Internet and similar, must be kept at a distance from BIOSTIM SPT, to avoid the risk of interference. The recommended minimum separation distance depends on the output power of the RF device and the transmission frequency. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications systems and BIOSTIM SPT, using reference the table below.

Test Frequency (MHz)	Band (MHz)	Type of service	Modulation	Maximum power (W)	Distance (m)	Immunity Test Level (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0,3	28
710						
745	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
780						
810						
870	800 – 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
930						
1 720						
1 845	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
1 970						
2 450	2400 – 2570	Bluetooth, W LAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
5 240						
5 500	5100 – 5800	W LAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
5 785						

The immunity levels shown in the table are respected as long as the device is maintained at a distance of at least 30 cm from any possible source of RF disturbance.

BIOSTIM SPT, complies with all test levels with a distance >= 30 cm.

Note: BIOSTIM SPT, at the frequency of 450 Mhz is immune to the disturbance up to a level of 20V / m - To obtain the required immunity level (28 V / m) the distance from the source of disturbance must be increased to not less than 42cm.

## 8. SYMBOLS

	<b>Class II appliance:</b> Appliance in which protection against electric shock does not rely on basic insulation only, but includes additional safety precautions such as double insulation.
	<b>Appliance with BF type applied part:</b> appliance with a specific degree of protection against electrical hazards, specifically regards admissible leakage current, and with an F type (floating) applied part isolated from the rest of the appliance (coil).
	<b>Continuous current:</b> on the generator data plate indicates the required power supply, on the power supply data plate indicates the provided power supply.
	<b>Alternate Current (AC):</b> on the power supply data plate indicates the required power supply.

	<b>Attention, consult the enclosed documents:</b> symbol informing the user to consult the documents supplied with the appliance, as the user manual, for correct understanding and/or use of the part marked with this symbol.
	<b>Operating Instructions:</b> refer to the instructions supplied with the device for better use.
	<b>Manufacturer's name and address:</b> manufacturer's name and address.
	<b>Manufacturing date:</b> manufacturing year of the device.
	<b>Keep dry:</b> the device is not protected against the penetration of liquids and must be stored and used in dry.
	<b>Model:</b> This symbol, specified on generator's and battery's data plates, indicates the generator's and battery's model.
	<b>Serial Number:</b> This symbol, specified on the generator's data plate, indicates the generator's serial number.
	<b>Batch code:</b> This symbol indicates the battery batch code.
	<b>Separated collection:</b> IGEA devices and their applied parts cannot be disposed of with urban waste but must be disposed of separately according to requirements specified by local authorities.
	<b>Follow the instructions for use:</b> This symbol indicates that you must read the instructions before you start using the device.
	<b>CE symbol:</b> on the device's data plate indicates that the device conforms to European Directive for medical devices 93/42/EEC and its revised version. The CE symbol is followed by the identification number of the notified body that checks conformity to the Essential Requirements of the Directive for applying the mark on the product.
	<b>Limits of atmospheric pressure for the transport and storage:</b> symbol on the outer packaging of the system.
	<b>Limits of relative humidity for the transport and storage:</b> symbol shown on the external packaging of the system.
	<b>Temperature limits for the transport and storage:</b> the symbol found on the outer packaging of the system.
	<b>Fragile handle with care:</b> symbol on the outer packaging of the system.

Note: Any symbol not described in this manual is to be considered for exclusive use of the manufacturer to facilitate the internal management of the products, but is not intended to supply information to the user.

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