




Adjustable Electrodes

Over the last 5 years over 5,000 patients have been successfully treated with **the current CLINIPORATOR disposable electrode**. Over this time we've listened to our customers and this highlighted a number of development areas. In conjunction with surgeons and utilising the latest technology we have developed this fully adjustable, disposable electrode (patent pending), with the aim to ensuring ease of use and safety for the user and patient.



- › Better control during skin penetration
- › Avoids convergence of the needles
- › Depth of needle penetration can be pre-set



- › Totally disposable
- › Ready for use
- › No recapping required
- › Longer connection lead



Adjustable Electrodes



Adjustment in 5 mm increments 5-40 mm, 1 electrode treats different lesions without change



BOX CODE (5 PCS)	DESCRIPTION
IGOE070	adjustable electrode, linear configuration (Linear) length 30 mm (30), non-insulated needles (Standard) Double packaging Bag + Blister
IGOE080	Adjustable electrode, linear configuration (Linear) length 40 mm (40), partially insulated needles (insulated) Double packaging Bag + Blister



Adjustment in 5 mm increments 5-40 mm, 1 electrode treats different lesions without change



BOX CODE (5 PCS)	DESCRIPTION
IGOE050	adjustable electrode hexagonal configuration (Hexagonal) length 30 mm (30), non-insulated needles (Standard) Double packaging Bag + Blister
IGOE060	adjustable electrode hexagonal configuration (Hexagonal) length 40 mm (40), partially insulated needles (insulated) Double packaging Bag + Blister